

Digital Tele Thermometer

461

Features

- Extremely useful for pharmaceutical Labs
- High accurate and easy to operate
- 3½ Digit Display
- 6 Temperature Probes



1 Year
Warranty

Digital Telethermometer Model 461 is an ideal, low cost, accurate and easy to operate instrument useful for Pyrogen testing in pharmaceutical industries. The instrument has 6 channels so that temperature of 6 different rabbits can be monitored simultaneously. Special temperature probes are provided which can be easily inserted in the anus of the rabbits. Any of the 6 channels can be selected by a rotary switch provided on the front panel.

The instrument is extremely useful for pyrogen testing



Digital μ le Thermometer

Technical Data

Range	0 to 60°C
Resolution	0.1°C
Accuracy	0.2% FS \pm 1 Digit
Input	RTD (PT-100) or IC Sensor
Channels	6 channels, switch selectable
Repeatability	\pm 1°C
Display	3 1/2 Seven Segment bright red led
Power Supply	230V \pm 10%AC, 50 Hz
Dimensions	Inst: 275 X 175 x 76 mm Box: 325 x 300 x 230 mm
Weight	2.5 Kg (Approx.)
Accessories	Pyrogen Probes (Set of 6 Probes), Operation Manual, Dust Cover

Specifications are subject to change due to continuous improvements

Product Range

Semi Auto Analyser, Single Beam and double Beam UV/VIS Spectrophotometer, Flame Photometers, Haemoglobin Meter, Digital pH, Meter, Digital Photo Fluorometer, Digital Conductivity Meter, Digital Photo Colorimeter, Turbidity Meter, Auto KF Titrimer, Digital DO Meter, Digital Conductivity Meter, Digital TDS Meter, Temperature Indicator, Digital Telethermometer, Water & Soil Analysis Kit, Digital Salinity Meter, Digital Colony Counter, Portable pH Meter, Portable Conductivity Meter, Portable TDS Meter, Portable DO Meter, Digital Tablet Disintegration Apparatus, Digital Friability Test Apparatus Tablet Dissolution Test Apparatus, Melting Point Apparatus, Student Microscope, Binocular microscope, Medical Microscope, Co-Axial Microscope.