

Dissolved Oxygen Oxygen Meter Meter

831



Features

- Highly Stable & Accurate
- D.O. & Temperature Measurements
- 3½ Digit LED Display
- Gold Silver Amperometric Probe

1 Year
Warranty

Digital D.O. Meters, model 831 are available in bench and portable models. These are extremely useful for agriculture and soil analysis laboratories, fisheries, water quality control in boiler feed water, water works department, breweries, water purification plants, chemical & pharmaceutical industries etc.



Dissolved Oxygen Meter

Technical Data

MODELS	831
PRODUCT	Portable Digital DO Meter
DISPLAY	3½ Digital LCD
RANGE DO TEMPERATURE	0 to 20 ppm ---
RESOLUTION DO TEMPERATURE	0.1 ppm ---
CALIBRATION	Manual
ACCURACY DO TEMPERATURE	± 0.2 ppm ---
TEMPERATURE COMPENSATION	Manual 0 to 50°C
STIRRER	—
DIMENSIONS	180 X 85 X 30mm
WEIGHT	250gm
POWER	9 V DC Battery
ACCESSORIES	Gold/Silver Amperoric DO Probe, Set of DO Membranes, Dust Cover, Operation Manual, Temperature Probe (801), Magnetic Stirrer & Capsule (801)

Product Range

Semi Auto Analyser, Single Beam and double Beam UV/VIS Spectrophotometer, Flame Photometers, Haemoglobin Meter, Digital pH, Meter, Digital Photo Fluorometer, Digital Conductivity Meter, Digital Photo Colorimeter, Turbidity Meter, Auto KF Titrimeter, Digital DO Meter, Digital Conductivity Meter, Digital TDS Meter, Temperature Indicator, Digital Telethermometer, Water & Soil Analysis Kit, Digital Salinity Meter, Digital Colony Counter, Portable pH Meter, Portable Conductivity Meter, Portable TDS Meter, Portable DO Meter, Digital Tablet Disintegration Apparatus, Digital Friability Test Apparatus Tablet Dissolution Test Apparatus, Melting Point Apparatus, Student Microscope, Binocular microscope, Medical Microscope, Co-Axial Microscope.

Note:- Packaging in shock proof Styrofoam packing inside an attractive and thick plywood storage box