

Dual Channel Flame Photometer

391/392

Advantages

- Mixing Chamber Designed For Stable Flame
- Simple & Easy To Operate
- Narrow Band Interference Filters
- Accurate Na, K, Ca & Li Determination
- Dual Channel For Simultaneous Observation Of Na & K With Single Aspiration (391)
- Concentric Non-Corrosive Nebulizer
- Auto Ignition



1 Year
Warranty

Digital Clinical Flame Photometer Model 391 is a Dual Channel Flame Photometer used for simultaneous observation of Na & K in blood serum (Plasma), Urine & other type of biological fluids.

The instruments are supplied in two units viz. Main Unit and Compressor Unit. The Main Unit comprises of digital display, detector, filters and nebulizer. The compressor unit supplies the compressed air to the Main Unit.

Air at a given pressure is passed into the atomiser. The suction produced by it draws the solution of the sample into the nebulizer and is sprayed as fine mist into the flame. The light emitted by the element at the characteristic wavelength is isolated by the narrow band interference filters, which then falls on to the silicon photodiode. After signal conditioning of the photodiode output, it is sent to digital converter and the results are displayed directly in ppm/meq value. The narrow band interference filters are mounted on a rotating disc fitted within the Main Unit. Appropriate filter can be selected by rotating this disc with the thumb. The oil free compressor with regulator maintains a stable airflow.



Dual Channel Flame Photometer

Technical Data

COMPRESSOR UNIT

Air Supply : By oil free mini compressor unit with pressure regulator with load and over heating protection

Combustion Gas LPG : Controlled by precision regulator (Domestic Gas Cylinder)

Power : 230 V \pm 10% AC, 50 Hz

Dimensions(LxBxH) : 290 x 255 x 210 mm (Approx.)

Weight : 8 Kg. (Approx.)

ACCESSORIES

Operation & Service Manual - 2 Nos

Ignition Filament - 1 No

Capillary Tube Set - 1 Set

Cleaning Wire Set - 1 Set

Dust Protection Cover - 1 No

Nose Plier - 1 No

Technical Data		
MODEL	391	392
RANGE	Serum Na : 0-200 meq/l K : 0-10 meq/l Ca : - Li : 0-2 meq/l (Optional)	Urine 0-200 meq/l 0-100 meq/l 0-10 meq/l (Optional) - 0-100 ppm 0-100 ppm 0-100 ppm
READOUT	±2% upto 100 meq ±5% above 100 meq	2% up to 100 ppm 5% above 100 ppm
ACCURACY	Dual 2½ digit bright red LED display	
IGNITION SYSTEM	In-built electronic Ignition by press of a switch	
DIMENSIONS (LxBxH)	375 x 245 x 220 mm (Approx.)	
MEASUREMENT SYSTEM	2 elements on dual LED display	
REPEATABILITY	± 2 Counts	
DETECTOR	Silicon Photodiode	
FILTERS	Narrow band interference glass filters (Na & K)	
AUTOMISER	S.S. non-corrosive & easily accessible, Black bakelite axial flow type	
FLAME SYSTEM	LPG & dry oil free air	
WARM UP TIME	10 minutes	
POWER	230 V±10% AC, 50 Hz	
WEIGHT	Main Unit : 7.5 Kg. (Approx.) Compressor : 8 Kg. (Approx.)	

Rugged & accurate instrument to determine concentration of Na, K, Ca & Li in any solution.

Product Range

Double Beam UV/VIS Spectrophotometer, Single Beam UV/VIS Spectrophotometer, Flame Photometer, Haemoglobin Meter, Photo Colorimeter, pH Meter, Conductivity Meter, TDS Meter, DO Meter, Salinity Meter, Turbidity Meter, Nephelometer, Colony Counter, Karl Fischer Moisture Titrator, Fluorometer, Dissolution Test Apparatus, Disintegration Test Apparatus, Friability Apparatus, Melting Point Apparatus, Tap/Bulk Density Apparatus, Hardness Tester, Leak Test Apparatus, Flocculator/Jar Test Apparatus, Telethermometer, Temperature Indicator, Water & Soil Analysis Kit, Portable pH Meter, Portable Conductivity Meter, Portable TDS Meter, Portable DO Meter.

Other Flamephotometers available : Microprocessor Flame Photometer Model 381,391, 1381 & 1382

**ENVIRONMENTAL
AND SCIENTIFIC
INSTRUMENTS CO.**
AN ISO 9001 : 2015 CERTIFIED CO.

